

There are several criteria that are important to the professional chemist and the laboratory professional. These include sample preparation, reproducibility, high sensitivity, and versatility. Infrared (IR) analysis coupled Attenuated total reflectance (ATR) accessories is an excellent way to achieve these goals. This application describes how FTIR with ATR was used to analyze the polymers found in an automobile.

Introduction

Polymers are present in daily life. From the wrappings around the food we eat to the sacks used for the garbage, polymers are practically everywhere. The automobile is no exception. Polymers comprise greater than 40% of a modern automobile; the carpet, seats, padding, knobs, switches, dashboard and cosmetic coverings are composed of one polymer or another. This paper will describe the analysis of the various polymers in an automobile immediately accessible to investigation.




Figure 1. Transmission and ATR spectra of polystyrene.

The traditional method of analysis for polymers is to make a thin film of the polymer and collect an infrared transmittance spectrum. Sample preparation can be time consuming, film thickness inconsistent and the use of solvents hazardous. Also, not all polymers can be dissolved in a solvent, heated or flattened enough to make an adequate thin film for analysis with infrared spectroscopy. ATR accessories require no sample preparation and greatly simplify the collection of infrared (FT-IR) spectra.

The solid sample is placed onto the ATR crystal, pressure is applied, and the sample spectrum is collected. The sample is removed from the crystal surface and the accessory is ready to collect additional spectra. Faster and less complicated than using prepared thin films, the ATR method allows recovery of the original sample for other analysis methods.




Figure 2. Spectrum from vehicle carpet fibers

The resulting sample data can also be searched against a digital database of ATR spectra for positive identification. Despite changes in relative peak intensity of the absorption bands, due to the internal mechanism of ATR accessories¹, spectra can also be compared to transmission data. As an example, Figure 1 is a plot of the transmission and ATR spectra of polystyrene.

Experimental Results

Spectra were collected using a Jasco FTIR-460 Fourier transform infrared spectrometer with a diamond micro-ATR accessory. Sample sizes of 1-4 mm² were collected from various parts of the vehicle. The individual polymer samples were placed onto the ATR crystal and the anvil was hand-tightened to apply pressure. No sample preparation was necessary to collect the spectra.




Figure 3. Fibers from the seat covering and fabric of the vehicle roof.

FT-IR spectra of 64 scans at 4 cm^{-1} resolution were co added and averaged to obtain the single-beam background and sample spectra. To provide identification of the samples, the spectra were searched versus several library collections of polymer spectra.




Figure 4. Representative spectra of the Polymer samples collected from various molded side panels and trim.

Figure 2 is a spectrum of fibers comprising the nap of the vehicle carpet and is identified as Nylon-6. Fibers from the seat covering and fabric of the vehicle roof were analyzed and identified as Poly(ethylene terephthalate). Figure 3 is a representative spectrum of these fibers.




Figure 5. Spectrum of the trim from the radio console

Polymer samples collected from the various molded side panels and trim were identified as the complex polymer Poly(ethylene:propylene:diene). One of four sample spectra can be found as Figure 4.




Figure 6. Spectrum of the red polymer found in the taillight.

Figure 5 is a representative spectrum of samples collected from the trim around the radio console and a portion of the taillight lens, subsequently identified as Bisphenol-A-polycarbonate. The red polymer of the taillight was identified as Poly(methyl methacrylate) and is plotted as Figure 6.




Figure 7. Dashboard knobs identified as (80% styrene) Poly(styrene:acrylonitrile)

Several samples collected from the vehicle knobs, switches and the radio console were identified as a copolymer of styrene and acrylonitrile. The peak at 2235 cm^{-1} is representative of the nitrile functional group in polyacrylonitrile.

Conclusions

FTIR spectra collection is greatly simplified using a diamond single reflection ATR. ATR is rapid, repeatable and reliable for the characterization of polymers in any configuration. The analysis method is non-destructive and can be used to collect data from a minimal amount of sample.

Acknowledgements:

Specac Inc, 500 Technology Court, Smyrna, GA 30082
Tel: 800 447 2558 www.specac.com, for the ATR.

Contact:
Jasco, Inc.
8649 Commerce Dr. Easton, MD 21601
800-333-5272 Fax: 410-822-7526
www.jascoinc.com