

## System Evaluation of Spectrofluorometer FP-8500 with Optical Fiber

### Introduction

The OBF-832 optical fiber interface is used to connect an optical fiber for remote sampling of fluorescence materials with an FP-8000 Series fluorescence spectrophotometer. By using an optical fiber, measurement can be made by holding the probe close to the sample. Optical fibers can be used for a variety of measurements, such as when the sample is larger than the sample chamber, to follow an in-situ reaction, or for measurement in hazardous environments, such as high/low temperature and high pressure.



**FP-8500**  
Spectrofluorometer


In this application note, a comparison was made by measuring samples using an FP-8500 fluorescence spectrophotometer with both a FLH-809 film holder and an optical fiber with the OBF-832 optical fiber interface. To check for spectral accuracy the fiber probe was corrected for quantum efficiency using a calibrated light source with a validated spectral emission.

### Keywords

Spectrofluorometer, Optical fiber, Fluorescence spectrum, fluorescence spectrophotometer

### Measurement System and Condition

FP-8500	Fluorescence spectrophotometer
OBF-832	Optical fiber unit
Calibration	ESC-842 Calibration light source / standard white reference plate
Samples	Red and Yellow reference cards


**Figure 1.** Measurement system

Measurement Condition	
Measurement Mode	Fluorescence
Excitation Bandwidth	10 nm
Emission Bandwidth	10 nm
Excitation Wavelength	310 nm
Measurement Range	400 - 800 nm
Scan Speed	100 nm/min
Data Interval	0.5 nm
Response	1 sec
PMT Voltage	350 V (yellow card), 550 V (red card)


## Measurement Results

Figures 2 and 3 show the measurement results. The profiles of the spectra of each colored sample measured using both the optical fiber and FLH-809 film holder were consistent with each other.

These results demonstrate that the spectral accuracy obtained using both the fiber probe and film holder yield comparable data.



**Figure 2.** Normalized fluorescence spectra of the yellow card



**Figure 3.** Normalized fluorescence spectra of the red card