

## Determination of the Relative Quantum Yield of Rhodamine B

### Introduction

Different molecular and environmental conditions not only effect whether a molecule will fluoresce or not, but can also determine the intensity of the emitted fluorescence radiation. A molecule's efficiency to fluoresce is described by its quantum yield and is defined as the ratio of the number of photons absorbed to the number of photons emitted by the sample. There are two methods for measuring the fluorescence quantum yield: the absolute method and the relative method. The absolute method directly obtains the quantum yield by detecting all sample fluorescence through the use of an integrating sphere. The relative method compares the fluorescence intensity of a standard sample with the fluorescence intensity of an unknown sample to calculate the quantum yield of the unknown sample. Therefore, the obtained results depend on the accuracy of the standard sample's quantum yield value.



**FP-8500**  
Spectrofluorometer

This application note will demonstrate how to obtain the relative quantum yield of Rhodamine B using fluorescein as a standard sample.

### Keywords

FP-8500, FUV-803, Fluorescence, Quantum Yield, Relative Method

### Experimental

In order to calculate the relative quantum yield of an unknown sample, the following information is required: the quantum yield of the standard sample, the absorption spectra of the samples, and the area of the samples' emission spectra after spectral correction. Additionally, when the solvent of the standard sample is different from that of the unknown sample, the average refractive index value is needed. A dilution ratio is also required if the standard and/or unknown samples are diluted.

By obtaining the aforementioned parameters, the relative quantum yield of unknown sample,  $\Phi_x$ , can be calculated using the following equation:

$$\Phi_x = \Phi_{st} \cdot \frac{A_{st} \cdot F_x \cdot n_x^2 D_x}{A_x \cdot F_{st} \cdot n_{st}^2 D_{st}}$$

where  $\Phi_{st}$  is the quantum yield of the standard sample,  $A_{st}$  and  $A_x$  are the absorbance of the standard and unknown samples, respectively,  $F_{st}$  and  $F_x$  are the areas of the standard and unknown emission spectra, respectively,  $n_{st}$  and  $n_x$  are the refractive indices of the standard and unknown samples, respectively, and  $D_x/D_{st}$  is the dilution ratio.

Prior to obtaining the relative quantum yield calculations, rhodamine B was used to correct the excitation spectrum and a calibrated halogen light source was measured from 450 to 700 nm to correct the emission spectrum.

Measurement Conditions		
	Absorbance	Emission
Excitation Wavelength		500 nm
Excitation Bandwidth	2.5 nm	5 nm
Emission Bandwidth	10 nm	5 nm
Scanning Speed	200 nm/min	200 nm/min
Data Interval	1 nm	0.5 nm
Response Time	0.5 sec	0.5 sec
PMT Voltage	230 V	430 V

For absorbance measurements, 5400  $\mu\text{g/L}$  and 72  $\mu\text{g/L}$  of fluorescein and rhodamine B were dissolved in ethanol, respectively. For fluorescence measurements, 7200  $\mu\text{g/L}$  and 36  $\mu\text{g/L}$  of fluorescein and rhodamine B were dissolved in ethanol, respectively.

## Results

The absorption spectra of fluorescein and rhodamine B were measured using the [Absorbance Measurement] program and FUV-803 Absorbance measurement cell block. The results are illustrated in Figure 1 and indicate that the absorbance of fluorescein at 450 nm was 0.490 and the absorbance of rhodamine B at 500 nm was 0.225.




Figure 1. Absorption spectra of fluorescein (blue) and rhodamine B (purple).




Figure 2. Fluorescence emission spectra of fluorescein (blue) and rhodamine B (purple).

The emission spectra of rhodamine B and fluorescein were then measured and are shown in Figure 2. In order to prevent inner filter effects or reabsorption of the emitted fluorescence radiation, the sample solutions were diluted so that the sample absorbance was less than 0.02.

The area under fluorescein's emission spectrum was 727204 and the area under rhodamine B's emission spectrum was determined to be 243513. Since ethanol is used as the solvent for both the standard and unknown solutions, the average refractive index of the solvent is not required. However, when the solvents used are different, the published value is necessary. The standard sample solution was diluted by 75 times and the unknown sample 100 times, which makes the dilution ratio 1.33. The relative quantum yield of Rhodamine B was then calculated by applying the obtained parameters to the equation above. The calculated relative quantum yield of rhodamine B is 92% and within range of the published values of 69 - 97%<sup>1</sup>.

## Conclusion

This application note demonstrates that the FP-8500 can be used to easily determine a sample's relative quantum efficiency from a standard sample.

## References

1. Kazuhiko Kinoshita and Koshin Mihashi. Fluorescence measurements - Applications for Biochemical Sciences. The Spectroscopic Society of Japan, Measurement Method Series 3, Japan Scientific Societies Press, 1983.