

Time Course Measurement with FTIR

Introduction


FTIR offer three methods for time course measurement: interval measurement, rapid-scan, and step-scan; the time resolution is 1 sec, 50 msec, 5 msec (optional: 10 nsec) respectively.

Here we will review the basic features and measurement for each method. Each method uses a different type of scanning method for the moving mirror in the interferometer. “Interval measurement” is performed with regular scanning and “rapid-scan measurement” uses a faster scan rate of the moving mirror. Using “step-scan measurement,” the moving mirror stops at each data-sampling point and perturbed at each of the points.




FT/IR-4000 and 6000 Series Spectrometers

Measurement Conditions		
Interval Measurement¹	Time Res.: 1 sec ~	Application: structural transition of a sample, monitoring of gas concentration
Measurement Ex.: Measure the transition of CO ₂ concentration in room air. (Fig. 1) Interval measurement is effective for long-time monitoring of gas analysis. ¹ JASCO interval measurement program can be used for up to 60001 scans or 24 days.		
Rapid-Scan Measurement	Time Res.: 50 msec ~	Application: photo-polymerization reaction, orientation relaxation of polymer film
Measurement Ex.: The curing process of a UV cured resin was measured, a photo-polymerization reaction. (Fig. 2) The peak intensity transition at 1637 cm ⁻¹ (bottom of Fig. 2) is attributed to the C=C stretch. The peak intensity rapidly decreased responding to the UV irradiation and continued to decrease moderately indicating that the initial reaction and reaction process for photo-polymerization was observed very accurately.		
Step-Scan Measurement²	Time Res.: 5 msec ~ ³	Application: relaxation of liquid crystal orientation
Measurement Ex.: The relaxation process of the orientation of a liquid crystal in response to an electric field was measured. (Fig. 3) The peak intensity transition at 2925 cm ⁻¹ (bottom of Fig. 3), is attributed to the C-H stretch. The peak intensity fluctuated responding to the power voltage being applied and stopped, and decreased in two phases - rapidly and then moderately. It is known that the liquid crystal in the area near the electrode has a faster orientation relaxation than the bulk and this difference can be observed.		
² limited to the materials with repetitive response. ³ 10 nsec ~ is optionally available.		




3D Spectrum




Peak intensity transition at 2360 cm⁻¹

Figure 1. CO₂ concentration transition




3D Spectrum




Peak intensity transition at 1637 cm⁻¹

Figure 2. Curing process of a UV cured resin



3D Spectrum



Peak intensity transition at 2925 cm⁻¹

Figure 3. Relaxation process of liquid crystal orientation